

# Histone Maleimide Labeling

## Materials:

- Atto488 maleimide (Sigma, 28562)
- 5ml HiTrap desalting column (GE, 29-0486-84)
- TCEP
- Lyophilized histone

## Labeling Buffer:

	50 mL	Final Concentration
1M Tris, pH7.0	1mL	20mM
GuCl (MW: 95.5g/mol)	33.4g	7M
0.25M EDTA	1mL	5mM

## Procedure:

1. Bring up TCEP to 0.5M, adjust pH to 7.5.
2. Bring up histone in labeling buffer to 0.5mM final concentration.
  - for H2A 120C, 3.6 mg of histone goes into 516uL labeling buffer
  - epsilon = 0.317, MW = 14081g/mol
1. Add TCEP to histone to final concentration of 0.7mM.
2. Incubate at RT for 2 hrs.
3. Dissolve Atto488 maleimide in DMSO to 100mM, vortex and spin down to completely dissolve. Wrap tube in foil and store at -20C.
  - for 1mg Atto488, dissolve in 9.37uL DMSO
  - Atto488 maleimide MW: 1067g/mol
1. Add Atto488 maleimide to histone:
  - for 5ml desalting column, optimal load volume is 100uL-1.5mL; 200uL reaction volumes work well.

- alternative desalting method is to dialyze labeling reaction into labeling buffer in which case you don't need to worry about reaction volumes; dialysis desalting method seems to work well for Atto488 and Atto565

	50 mL	Final Concentration
100mM Atto488	1.2uL	0.6mM
0.5mM histone	40uL	0.1mM
Water	158.8uL	

Incubate 3 hr at RT, then O/N at 4C

1. Equilibrate HiTrap desalting column with labeling buffer at 4C.
  - a. Load sample at 1ml/min, collect 200uL fractions.
  - b. Observe 2 peaks, 1<sup>st</sup> narrow and 2<sup>nd</sup> broad
2. Run SDS-PAGE and visualize with Typhoon imager in Atto488 channel, then coomassie stain.
  - a. 1st peak corresponds to protein (labeled and unlabeled elute together)
  - b. 2nd peak corresponds to free dye
3. Assess labeling efficiency by NanoDrop (delta between molarity calculated from Atto488 absorbance and molarity calculated from histone absorbance)