

# Preparing Lipids

Components, from Avanti Polar Lipids:

**DOPC:** 1,2-dioleoyl-sn-glycero-3-phosphocholine, as 1g powder

**PEGylated lipids:** 18:1 PEG2000 PE

1,2-dioleoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-2000] (ammonium salt), as 100mg powder

**biotinylated lipids:** 18:1 Biotinyl Cap PE

1,2-dioleoyl-sn-glycero-3-phosphoethanolamine-N-(cap biotinyl) (sodium salt), as 25mg powder

**Prepare the Master Mix (MM)** in the DOPC glass bottle. The final proportions are:

100mg/ml DOPC

10mg/ml PEG-2000 DOPE

.5mg/ml biotinylated DOPE

Use up the entire 1000mg of the DOPC, so the final volume of chloroform will be 10ml. This translates to 100mg total PEG-2000 DOPE and 5mg biotinylated DOPE.

- 1) Take out all components from the freezer and leave at room temperature for 30 minutes to 1 hour. This will prevent water condensation when the vials are opened.
- 2) Thoroughly clean a 10ml glass pipette.
- 3) Open the 1g DOPC bottle and add 100mg PEG-2000 DOPE directly into it. The PEG-2000 DOPE is, conveniently, available as a 100mg vial.
- 4) Then measure out and add ~5mg biotinylated PEG.
- 5) Rinse the glass pipette with chloroform. Dissolve the dry lipid mix in 10ml chloroform. Mix gently by swirling until all the dry lipid is dissolved.
- 6) Seal the bottle with plenty of parafilm and store at -20.

- 1) **To prepare lipids for experiments**, rinse a small glass vial with water then alcohol; dry it with nitrogen, and place it in the vacuum oven for 30 minutes or so.
- 2) Take the MM out of the freezer and let it come to room temperature to prevent condensation on opening.
- 3) Rinse a 1ml (or similar) Hamilton syringe (with a metal plunger!) with chloroform. Then measure out 200ul of the MM and transfer it to the cleaned glass vial.
- 4) Slowly evaporate the chloroform with a slow nitrogen stream. Try to get the lipids to dry in a film across the side of the vial, up to a height of a few cm. Make sure all the chloroform, including pockets caught under the lipid film, is evaporated. You can up the nitrogen pressure as the chloroform evaporates.
- 5) Place the vial under vacuum overnight. It can stay like this for several days if necessary.
- 6) Resuspend the lipids in 2ml lipid buffer at room temperature for at least an hour. It can go overnight.
- 7) Vortex and sonicate. Ask someone about sonication procedures for lipids—everyone has a different approach.